

**QUICKIE**<sup>®</sup>

**Pulse**<sup>™</sup>

Comfort and performance  
in a proven design



# QUICKIE® Pulse

## AVAILABLE IN 6 MODELS:

### Quickie Pulse 5

BC (K0848): 5 MPH, 2-Pole motor,  
No Power Option base with ASAP II Rehab Seat

CC (K0849): 5 MPH, 2-Pole motor,  
No Power Option base with Captain's Seat

### Quickie Pulse 6

BC (K0848): 6 MPH, 4-Pole motor,  
No Power Option base with ASAP II Rehab Seat

CC (K0849): 6 MPH, 4-Pole motor,  
No Power Option base with Captain's Seat

SC (K0856): 6 MPH, 4-Pole motor,  
Single Power Option Tilt (SPOT) or Power Recline with  
ASAP II Rehab Seat or 3.7 Recline Seat

MPC (K0861): 6 MPH, 4-Pole motor, Power Tilt and Recline  
with 3.7 Recline Seat



## POWER BASE DESIGNED FOR PERFORMANCE

- Mid wheel drive
- Reliable and compact for indoor accessibility
- Flexible suspension system for optimal outdoor performance
- Black or gray tire color options
- Available with 15 shroud color options
- Reliable 2-pole and 4-pole motors simulated a mean time of 20 years to failure during High Cycle Fatigue test

## ELECTRONICS OPTIONS

### VR2 Non-Expandable Controls

- Industry standard non-expandable electronics
- Dependable, highly flexible, and simple to use
- 4-way toggle switch for operation of tilt and power legrests without the need for expandable controls

### R-net Expandable Controls

- Assignable buttons for convenient control of common functions
- Optional full color display
- 8 uniquely programmable profiles for easy transitions between environments



## MULTIPLE SEAT FRAME OPTIONS



### Captain's Seat (CC)

- Fully adjustable seat size, legrest placement, and armrest width
- Convertible seat pan for flexibility in cushion selection
- Comfortable style
- Semi-recline and fold-down backrest
- Optional O2, crutch, and walker holders



### ASAP II Rehab Seat (BC & SC)

- In-field width, depth, armrest, legrest, and back height adjustability
- Manual reclining and fold-down backrest
- Available with back upholstery, SPOT Curved backrest, JAY® & ART SPO customizable backrests, and the JAY J3 Back
- Extensive seating and positioning capabilities
- Optional O2, crutch, and walker holders

### 3.7 Recline Rehab Seat (SC & MPC)

- Option of Recline Curved back or Integrated JAY J3 Back for seamless functionality with power recline and shear reduction system
- Low profile electronics enclosure shell for a clean look and compatibility with accessories
- Seat-to-floor heights as low as 16.3"
- Cantilever Reclining Armrests stay parallel to ground when reclined

## POWER SEATING OPTIONS

- NEW 3.7 Power Recline with Power Shear Reduction
- 50° of Power Tilt to optimize positioning
- 9" Power Seat Elevate for increased access to cabinets and eye-to-eye interaction with peers
- Power Elevating Legrests with Articulation
  - Swing-Away
  - Centermount with Single/Dual Footplate for angle adjustability as a single footplate or separate footplates



## Quickie Pulse Specifications

SPECIFICATIONS	PULSE 5 BC	PULSE 6 BC	PULSE 5 CC	PULSE 6 CC	PULSE 6 SC	PULSE 6 MPC
Drive Wheel Position	Mid	Mid	Mid	Mid	Mid	Mid
HCPCS Code	K0848	K0848	K0849	K0849	K0856	K0861
Speed	5.0 MPH	6.0 MPH	5.0 MPH	6.0 MPH	6.0 MPH	6.0 MPH
Overall Weight	235 lbs.	235 lbs.	250 lbs.	250 lbs.	260 lbs.	345 lbs.
Overall Width	24"	24"	24"	24"	24"	24"
Overall Length	34"	34"	34"	34"	34"	34"
User Weight Capacity	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.
Seat Widths	12" - 22"	12" to 20"	16", 18", 20"	16", 18", 20"	12" to 22"	14" to 22"
Seat Depths	12" - 22"	12" to 22"	15" to 22"	15" to 22"	12" to 22"	16" to 20"
Seat To Floor Height	16" - 20"	16" to 20"	16" to 18" (from pan)	16" to 18" (from pan)	16.3" to 21.6"	16.3" to 18.3"
Seating Frame Options	ASAP II Rehab Seat	ASAP II Rehab Seat	Captain's Seat	Captain's Seat	ASAP II Rehab Seat or 3.7 Recline	3.7 Recline
Transit Approved	ISO	ISO	WC-19*	WC-19*	Tilt: ISO/Recline: WC-19*	WC-19*
Warranty	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime

\* Conforms with ANSI/RESNA WC-19 standards. Please refer to product Owner's Manual for additional information. Owner's Manuals are available for download at [www.sunrisemedical.com](http://www.sunrisemedical.com)



## COMPLEMENT YOUR QUICKIE

Order your JAY® Cushion or Back shipped complete on-chair with any Quickie Pulse chair. Just check the boxes on the order form to save time and money!

### JAY® Back Options

Available on BC, SC and MPC Models



**J3™ Posterior**  
E2613



**J3™ Posterior & Lateral**  
E2615



**J3™ Posterior Deep Lateral/Contour**  
E2620



**Quickie Tension Adjustable Upholstery**



**JAY or ART SPO Backrest**  
E2613, E2617



**JAY SPOT Curved or Curved Recline Back**

### JAY® Cushion Options



**Ion™**  
E2607, E260



**Easy™**  
E2607, E2608



**Union™**  
E2307, E2608



**X2™**  
E2622



**J2™/J2™ Deep Contour**  
E2622, E2623,  
E2624, E2625



**Fusion™**  
E2622, E2623



**J3™**  
E2622, E2623,  
E2624, E2625



Customer Service: 800.333.4000  
[www.SunriseMedical.com](http://www.SunriseMedical.com)

©03.2014 Sunrise Medical (US) LLC  
110300 Rev. G