

TheraPedic™ Positioning Restraint System

THE **FIRST CHOICE**
OF PARENTS AND
THERAPISTS FOR
ASSISTING IN THE
SAFE TRANSPORT OF
INDIVIDUALS WITH
PHYSICAL DISABILITIES

Available Colors



Includes
(See page 4 for details)


- Adjustable Head Pads
- Shoulder Harness Pads
- Crotch Strap Pad
- Lateral Positioning Pads
- 5-Point Safety Harness
- Tether Strap Kit
- Locking Clip for Tether
- Easy to Follow Installation Instructions



- Features**
- Two sizes available, accommodating users weighing 20 - 130 lbs. and up to 5'6" tall
 - Durable, machine washable, comfortable padded cover
 - Vacuform ABS plastic shell and chromed steel side tubes
 - Standard 5-point harness, unlike many other seats that utilize the automobile's lap and shoulder restraint
 - Meets Fire-Retardant Code 302 for automotive upholstery
 - Latex and silicone free
 - Crash tested¹

- Benefits**
- Fits in most cars, vans and school buses
 - Fits most airplane seats and is certified for use in aircraft²
 - Lightest restraint system for use in vehicle, weighing much less than most comparable restraint systems
 - Multiple accessories allow restraint systems to grow with the child

Accessories (See page 4 and 5 for details)



MODEL NUMBER	2000	2500
Product Weight	11 lbs.	18 lbs.
Weight Capacity	20 - 102 lbs.	40 - 130 lbs. ³
Height Capacity	Up to 60	54 - 66
Shoulder Heights when Sitting	11½, 14½, 17½, 20½	13½, 17, 20½, 24
Depth to Seat Belt	7	9
(a) Seat to Top Height	29	33
(b) Seat Depth	10	12
Seat Depth with Extender ⁴	12½ or 15	16
(c) Inside Width at Front	15½	19½
Inside Width at Back	13	17½
Outside Width	19	23
Outside Depth	11	12
Seat to Back Angle	88°	85°

Note: Dimensions are in inches unless otherwise noted

¹ a. **MODEL #2000** meets or exceeds crash testing requirements in the US (FMVSS #213), Canada (NCSVA #213.3), and Europe (ECE R44.03)

b. **MODEL #2500** meets or exceeds crash testing requirements in the US (FMVSS #213) and Canada (NCSVA #213.3)

² **MODEL #2000** only

³ User must have a minimum of 13½" shoulder height to use this product

⁴ Seat depth increases with use of the Seat Depth Extender. See page 4 for details.

Accessories for TheraPedic™ Positioning Restraint System



ShurShape™ Custom Inserts

A unique way to provide custom molded seating for our TheraPedic Positioning Restraint System. Easily mold and capture a shape and then deliver the finished product to the user in one hour without a simulator.

MODEL #2030 - For use with **MODEL #2000**
MODEL #2530 - For use with **MODEL #2500**



Block Style Positioning Pads

Additional pads for special positioning needs. Available in four packs. For use with all TheraPedic Positioning Restraint Systems. Block style, two pads 2½" thick and two pads 1½" thick.

MODEL #2008 (Blue or Rose)



Swing Away Abductor

Improves hip alignment and stabilizes seating posture.
 Dimensions: 3½" W x 4" H x 3" L

MODEL #2015¹ - For use with **MODEL #2000**
MODEL #2018¹ - For use with **MODEL #2000** with 2½" extender
MODEL #2025
MODEL #2019¹ - For use with **MODEL #2000** with 5" extender
MODEL #2026
MODEL #2515¹ - For use with **MODEL #2500**
MODEL #2519¹ - For use with **MODEL #2500** with 4" extender
MODEL #2525



Lateral Positioning Pads

Four additional pads to improve lateral positioning. Inserted under the cover.

MODEL #2010 - For use with **MODEL #2000**
MODEL #2510 - For use with **MODEL #2500**



Seat Depth Extender

Increases seat depth.

For use with **MODEL #2000**
MODEL #2025¹ - Increase seat depth by 2½"
MODEL #2026¹ - Increase seat depth by 5"
 For use with **MODEL #2500**
MODEL #2525¹ - Increase seat depth by 4"



Extensor Thrust Wedge

Changes seat to back angle. Increases height by 2" or 3"

For use with **MODEL #2000**
MODEL #2050 - 2" front height to 1" rear
MODEL #2055 - 3" front height to 1" rear

For use with **MODEL #2500**
MODEL #2550 - 2" front height to 1" rear
MODEL #2555 - 3" front height to 1" rear



BuckleGuard

Prevents child from releasing the shoulder harness. For use with all TheraPedic Positioning Restraint Systems.

MODEL #2075

¹When ordering this item for an existing TheraPedic Positioning Restraint System, a new seat cover upholstery (**MODEL #2002** or **MODEL #2502**) must be purchased; please specify color Blue or Rose